

KeaExpert release notes

Note: items in **red** may require modifications to your own sequences. All personal sequences should be recompiled after each software update.

Version 2.02.15 (21-April-2026)

- Added experiments 1PulseAmplitudeSweep and 1PulseDurationSweep.
- Added a simulation mode for certain experiments – see ‘Select Kea’ in the file menu to choose this option.

Version 2.02.14 (21-March-2026)

- Ported experiments SSE and T1Sat from Mouse-Kea. Added T1Sat-ie which uses interexperiment time rather than repetition time.

Initial release Version 2.02.12 (13 February 2026)

- Port of 2.02.12 SpinsolveExpert to KeaExpert.
 - Removal of Spinsolve specific features.
 - Port of most important Kea-NMR and Kea-Mouse experiments.
 - Modification of preferences/common parameters interface to allow setting these values since for the Kea these may change more often than on the Spinsolve.
 - A new spectrometer parameter viewer/updater designed for the Kea.
 - Some experimental Mouse lift commands added